CardioMPO™ Control Kit

Catalog Number: 8601

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CardioMPO[™] Control Kit

INTENDED USE

The PrognostiX *CardioMPO*[™] Control Kit is intended for use with the *CardioMPO*[™] Reagent Kit as an assayed quality control sample to monitor and evaluate the precision and accuracy of the *CardioMPO*[™] Test.

SUMMARY AND EXPLANATION

The use of quality control materials is indicated as an objective assessment of the precision of methods and techniques in use and is an integral part of good laboratory practices. The *CardioMPO*[™] Control Kit contains three levels of MPO provided ready-to-use in frozen liquid form for extended shelf life.

REAGENTS AND MATERIALS

Materials supplied with the kit:

***CardioMPO*[™] Controls, 3 vials, MPO levels 1-3, 0.5 mL each**

Three concentrations of human MPO in human lithium heparin plasma. Refer to individual Control vial labels for actual concentration ranges of Controls.

Materials required but not supplied:

CardioMPO[™] Enzyme Immunoassay Reagent Kit Cat. Number 7601

WARNINGS AND PRECAUTIONS

For *in vitro* diagnostic use.

For professional use only.

Contains Human Source Material. All Controls should be treated as potentially biohazardous material. Human plasma matrix was tested for presence of antibodies to Human Immunodeficiency Virus Type 1 (HIV-1) and Type 2 (HIV-2), Hepatitis B Surface Antigen (HBsAg), and antibodies to Hepatitis C (HCV) and found to be negative/non-reactive. Because no test can offer complete assurance that HIV-1, HIV-2, Hepatitis B or C virus, or other infectious agents are absent, this material should be handled using good laboratory practice and Universal Precautions to avoid contact with skin and mucous membranes.^{1,2}

REAGENT STORAGE

Unopened vials must be stored frozen at -15 to -20°C. The Controls are provided ready to use.

STABILITY

Unopened frozen vials:

See Expiration Date on the Control label.

Thawed vials:

Controls are stable for 8 hours at room temperature (20-26°C) and for 14 days when stored at 2-8°C. Thawed vials that are in use or stored for longer than 1 hour should be recapped tightly to prevent evaporation.

INDICATIONS OF INSTABILITY OR DETERIORATION OF CONTROLS

Store all *CardiaMPO*™ Controls at the labeled temperature and do not use past their expiration dates.

Do not use any Controls that have indications of discoloration, cloudiness or precipitation. Do not use any Controls that show any signs of leakage.

ASSAY PROCEDURE

CardiaMPO™ Controls are treated in the same manner as patient samples and run in accordance with the procedure found in the package insert for the *CardiaMPO*™ Reagent Kit.

Before use, equilibrate Controls to room temperature (20 – 26°C) for a minimum of 30 minutes. Thoroughly mix the contents of the vial while avoiding foaming.

ASSIGNMENT OF CONCENTRATION RANGES

The concentration ranges were derived from replicate analysis and are specific to a particular lot of product. The testing was performed using a representative sampling of this lot of Controls. Individual laboratory means should fall within the acceptable range; however, laboratory means may vary from the listed values during the life of these Controls. Variations over time and between laboratories may be caused by differences in laboratory technique and instrumentation. It is recommended that each laboratory establish its own means and acceptable ranges and use those provided only as guides.

EXPECTED VALUES

Refer to individual Control vial labels for actual concentration ranges of each Control.

REFERENCES

1. HHS Publication (CDC) 93-8395, "Biosafety in Microbiological and Biomedical Laboratories," (1993).
2. CLSI (formerly NCCLS), Procedures for the Handling and Processing of Blood Specimens; Approved Guideline-Second Edition. Document H18-A2 [ISBN 1-56238-388-4]. CLSI, 9540 West Valley Road, Suite 1400, Wayne, PA 19087-1898, 1999.

WARRANTY

This product is warranted to perform as described in its labeling and in Prognostix, Inc.'s literature for this product when used in strict conformity with published directions regarding the use of this product. Prognostix, Inc. makes no other warranty, express or implied, and DISCLAIMS ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. UNDER NO CIRCUMSTANCES WHATSOEVER SHALL PROGNOSTIX, INC. BE LIABLE FOR ANY INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

Manufactured for:
PrognostiX, Inc.
10265 Carnegie Avenue, Cleveland, Ohio 44106
Tel: (866) 358-9828, (216) 445-7469 (PGNX)
Fax: (216) 444-5335
info@prognostix.com
www.prognostix.com

Authorized Representative:

Emergo Europe
Molenstraat 15
2513 BH, The Hague
The Netherlands
Tel: +31 (0)70 345 8570
Fax: +31 (0)70 346 7299

This product is covered by U.S. Patent Application 20020164662.
CardioMPO™ is a pending trademark of PrognostiX, Inc.

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CardioMPO™ è un marchio in corso di registrazione di PrognostiX, Inc.