



Adeno-Associated Virus serotype 2 (AAV2) ELISA

Catalog Number	PRATV
Name	AAV2 Titration ELISA
Assay format	Solid phase sandwich ELISA method.
Kit Control	3 vials, lyophilized
Assay Feature	Fast, sensitive & reproducible titration for AAV2 wt virions, AAV2 Gene Shuttles or assembled and intact empty AAV2 capsids
Sample Types	Cell Culture & purified virus preparation
Sample Volume	100 µL
Incubation	3 hours at 37°C, 10-15 min at room temperature
Sensitivity	5x10 ⁷ to 1x10 ⁹ particles/mL
Storage	2 - 8 °C
Configuration	96 Determinations, 12x8 removable strips
Clinical Fields	Gene Therapy Tool
FDA Status	For Research Use Only

This product is sold for laboratory research use or further manufacturing only and should not be used for human therapeutic or diagnostic applications. The information presented is believed to be accurate; however, said information and products are offered without warranty or guarantee since the ultimate conditions of use and the variability of the materials treated are beyond our control. Nothing disclosed herein is to be construed as a recommendation to use our products in violation of any patents. Under no circumstances shall ARP American Research Products, Inc. be liable for damages, whether consequential, compensatory, incidental or special, strict liability or negligence, breach of warranty or any other theory arising out of the use of the products available from ARP American Research Products, Inc. Nothing contained herein warrants that the use of the products will not infringe on the claims of any patents covering the product itself or the use thereof in combination with other products or in the operation of any process.