

**Product:** Pegylated Granulocyte Colony Stimulating Factor, Human Recombinant  
**Catalog #:** 11-PGCSF  
**Amount:** 5 µg

---

**DESCRIPTION:** Recombinant Human PEG-G-CSF produced in E.coli is a single, non-glycosylated, polypeptide chain containing 175 amino acids and having a molecular mass of 18,800 dalton. rHuPEG-G-CSF is manufactured by attaching a 20,000 dalton methoxypolyethylene glycol propionaldehyde (mPEG-ALD) to the N-terminal amino acid of G-CSF.  
G-CSF is purified by proprietary chromatographic techniques.

**SOURCE:** *E.Coli* Expression System

**PURITY:** Greater than 95.0% as determined by:  
(a) Analysis by RP-HPLC.  
(b) Anion-exchange FPLC.  
(c) Analysis by reducing and non-reducing SDS-PAGE Silver-Stained gel.

**ENDOTOXIN:** Less than 0.1 ng/µg (IEU/µg) of rHuPEG-G-CSF.

**DIMERS & AGREGATES:** Less than 1% as determined by silver-stained SDS-PAGE gel analysis

**FORM:** Purified, in 10mM sodium acetate buffer PH= 4 and 5% mannitol.

**PROTEIN CONTENT:** Protein quantitation was carried out by two independent methods:  
1. UV spectroscopy at 280 nm using the absorbency value of 0.815 as the extinction coefficient for a 0.1% (1mg/ml) solution. This value is calculated by the PC GENE computer analysis program of protein sequences (IntelliGenetics).  
2. Analysis by RP-HPLC, using a calibrated solution of PEG-G-CSF as a Reference Standard.

**CONCENTRATION:** 0.1 mg/ml

**BIOLOGICAL ACTIVITY:** rHuPEG-G-CSF is fully biologically active when compared to standard. The ED50, calculated by the dose-dependant proliferation of murine NFS-60 indicator cells (measured by 3H-thymidine uptake) is < 0.1 ng/ml, corresponding to a Specific Activity of 1 x 10<sup>8</sup> IU/mg.

**STORAGE:** -20°C (aliquot), avoid repeated freeze and thaw cycles

This product is sold for laboratory research use or further manufacturing only and should not be used for human therapeutic or diagnostic applications. The information presented is believed to be accurate; however, said information and products are offered without warranty or guarantee since the ultimate conditions of use and the variability of the materials treated are beyond our control. Nothing disclosed herein is to be construed as a recommendation to use our products in violation of any patents. Under no circumstances shall ARP American Research Products, Inc. be liable for damages, whether consequential, compensatory, incidental or special, strict liability or negligence, breach of warranty or any other theory arising out of the use of the products available from ARP American Research Products, Inc. Nothing contained herein warrants that the use of the products will not infringe on the claims of any patents covering the product itself or the use thereof in combination with other products or in the operation of any process.

489 Common Street, Belmont, MA 02478, USA; Email: info@arpl.com or staff@arpl.com  
Phone: 800-832-2611 or 617-489-1120 Fax: 888-832-2613 or 617-489-5120